


Foot Shell Ordering Guide Freedom Innovations

Please use the following guide to determine which Freedom Innovations cosmetic foot shells are available for each prosthetic foot.


- Most products can be provided with a Generation 3 foot shell with a Cosmetic Attachment Plate (CAP). The table at the bottom lists exclusions.
- Sandal toe foot shells are available for most products. The table at the bottom lists exclusions. The foot module must also be ordered in a sandal toe configuration (CAPs are not available for sandal toe products).



Cosmetic Foot Shell: Standard, with CAP



Cosmetic Foot Shell: Standard, without CAP



High-profile: Without CAP




Sandal Toe: Without CAP




Kinnex: Without CAP




Kinnex: With Sandal Toe, without CAP



Maverick: Without CAP



Maverick: With Sandal Toe, without CAP



Exclusion List:

Prosthesis	Foot Shell			
	Cosmetic Foot Shell: Standard, with CAP	Cosmetic Foot Shell: Standard, without CAP	High-Profile, without CAP	Sandal Toe, without CAP
Freestyle	-	✓	-	-
LP Symes	-	✓	-	-
Runway	-	✓	-	✓
Renegade AT	✓	✓	✓	-
Senator	✓	✓	✓	-
Defender	-	✓*	-	-

✓ = prosthesis available with this foot shell

*Generation 1 and 2 only.