

## Universal Abdominal Binder



### Features and Benefits

- Neoprene binder with touch and close fastening and bound edges
- Universal sizing to reduce stock holding
- Adjustable to allow for stoma care

### Indications for Use

- Post-op support/ventral hernia

### Ordering Information

Size	Front Depth (cm)	Back Depth (cm)	Product Code
Universal	22	22	NB

## Universal Deep Abdominal Binder



### Features and Benefits

- Deep abdominal binder with simple touch and close fastening
- Adjustable to allow for stoma care
- Universal sizing to reduce stock holding

### Indications for Use

- Post-op support/ventral hernia

### Ordering Information

Size	Front Depth (cm)	Back Depth (cm)	Product Code
Universal	22	30	AB

## Universal Temporary Hernia Support



### Features and Benefits

- A separate 150mm diameter, 18mm thick soft hernia pad can be located as required
- Supplied with 1x hernia pad, additional pads can be purchased separately

### Indications for Use

- Abdominal hernia
- General posture support
- Lumbar strains
- Sciatica
- Lumbosacral strains



### Ordering Information

Size	Front Depth (cm)	Product Code
Universal	23	1043
Universal	30	1044

Accessories	Product Code
Hernia Pad	1041/P