

## BlueROCKER® Short



### Features and Benefits

- Specifically developed to accommodate individuals with a shorter stature, slightly broader-shaped tibia, and wider foot at the metatarsal area
- Tibia padding should always be worn for patient comfort
- Both models have a heel height of 7mm in all sizes

### Indications for Use

- Foot drop in combination with weak plantar flexors, and in some cases limited knee control
- May also be used for bilateral problems and partial foot amputations, most proximal level is Chopart



### Ordering Information

| Size   | Shoe Size (UK) | Shoe Size (EU) | Footplate Length (cm) | Height (cm) | Product Code |           |
|--------|----------------|----------------|-----------------------|-------------|--------------|-----------|
|        |                |                |                       |             | Left         | Right     |
| Small  | 1 - 2.5        | 32 - 35        | 21.5                  | 29.5        | 289951011    | 289952011 |
| Medium | 2.5 - 5        | 35 - 38        | 23.5                  | 32          | 289951012    | 289952012 |
| Large  | 5 - 8          | 38 - 42        | 25.5                  | 34          | 289951013    | 289952013 |

## BlueROCKER® 2.0



### Features and Benefits

- The BlueROCKER® 2.0 is the updated version of BlueROCKER® Original
- Identical in shape to the ToeOFF® 2.0 but offers more rigid orthotic control
- Developed primarily for bilateral patients and those with more involved pathologies
- The extra stability will improve both balance and posture and give the user greater security, especially individuals with weak quadricep muscles
- Usually the preferred orthosis to be used in conjunction with a socket and toe filler for management of partial foot amputations
- SoftKIT™ 2.0 comprises padding shaped to precisely fit the orthosis and comes supplied with two tibia pads
- ComfortKIT™ 2.0 comprises thicker padding, designed for the more sensitive skin, fits the orthosis precisely

### Indications for Use

- Footdrop
- Weakness or impairment in multiple leg muscle groups
- Impaired balance
- Limb proprioception
- Partial foot amputation



### Ordering Information

| Size    | Footplate Length (cm) | Height (cm) | Product Code |           |
|---------|-----------------------|-------------|--------------|-----------|
|         |                       |             | Left         | Right     |
| X-Small | 21                    | 36          | 289401010    | 289402010 |
| Small   | 23                    | 38          | 289401011    | 289402011 |
| Medium  | 24.5                  | 40.5        | 289401012    | 289402012 |
| Large   | 27                    | 43          | 289401013    | 289402013 |
| X-Large | 28.5                  | 45          | 289401014    | 289402014 |

## BlueROCKER® 2½



### Features and Benefits

- The BlueROCKER® 2 ½ has the same dynamic functions and stability as a BlueROCKER® 2.0 with the exception that the 2½ has a lower heel height, 7mm, to adapt to modern shoes
- The toe lift is also slightly lower
- Spacious toe box
- Has all MikroFIX™ pre-applied to simplify the fitting of the padding
- Durable and strong straps
- The straps can be applied for either inside or outside closure on the orthosis
- Adaptable to the heel height of the shoe

### Indications for Use

- BlueROCKER® is designed to manage footdrop in conditions such as:
  - Stroke
  - Multiple Sclerosis
  - Post Polio Syndrome
  - Muscular Dystrophy
  - Spinal Cord Injuries
  - Traumatic Brain Injuries
  - Guillian-Barre Syndrome
  - Charcot-Marie-Tooth
  - Myelomeningocele
  - Neuropathy
  - Cerebral Palsy
- Supports gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amputations
- Partial foot amputations, most proximal level is Chopart

### Ordering Information

| Size    | Shoe Size (UK) | Shoe Size (EU) | Footplate Length (cm) | Height (cm) | Product Code |           |
|---------|----------------|----------------|-----------------------|-------------|--------------|-----------|
|         |                |                |                       |             | Left         | Right     |
| X-Small | 1 - 2.5        | 32 - 35        | 21.5                  | 36          | 289701010    | 289702010 |
| Small   | 2.5 - 5        | 35 - 38        | 23                    | 38          | 289701011    | 289702011 |
| Medium  | 5 - 8          | 38 - 42        | 24.5                  | 40.5        | 289701012    | 289702012 |
| Large   | 8 - 10.5       | 42 - 45        | 27                    | 43          | 289701013    | 289702013 |
| X-Large | >10.5          | >45            | 28.5                  | 45          | 289701014    | 289702014 |

## Comfort Night Splint



### Features and Benefits

- Adjustable flexion straps provide a gentle stretching of the plantar fascia and achilles tendon
- Soft elastane lined foam interior is soft, washable and breathable to eliminate heat build-up
- Elastane lined foam padded straps maintain comfort across the flexor tendons
- The orthosis helps to prevent pain and further injury by limiting contractions of the soft tissues at night
- Lightweight, low profile shell is sturdy and breathable
- Toe wedge to increase the plantar fascia stretch
- Ambidextrous fit - left/right

### Indications for Use

- Plantar fasciitis
- Achilles tendonitis
- Plantar flexion contractures

### Ordering Information

| Size   | Shoe Size (UK) | Product Code |
|--------|----------------|--------------|
| Small  | <6             | UFO/S/E      |
| Medium | 6 - 9          | UFO/M/E      |
| Large  | 9+             | UFO/L/E      |


[More Information](#)  
[Wear and Care Instructions](#)