

Basic Foam Collar



Features and Benefits

- Symmetrical design provides comfortable support
- Suitable for night use
- Touch and close fastening
- Made from medium density foam covered with stockinette

Indications for Use

- Spondylitis
- Spondylosis
- Rheumatoid arthritis
- Osteoarthritis

Note: Made to measure service available, please provide neck circumference and front depth.

Ordering Information

Size	Neck Circ. (cm)	Front Depth (cm)	Product Code
X-Small	33 - 37	7	BFC/XS
Small	38 - 42	8	BFC/S
Medium	43 - 46	9	BFC/M
Large	47 - 51	10	BFC/L
X-Large	52 - 56	11	BFC/XL

Need a different size to fit your patient?

Made to Measure service is available for this product!
To order, please provide neck circumference and front depth.



More Information
Wear and Care Instructions

Shaped Foam Collar



Features and Benefits

- Symmetrical design provides comfortable support, suitable for night use
- Increased shaping to reflect underlying anatomy and provide an improved fit
- Touch and close fastening
- Made from medium density foam covered with stockinette

Indications for Use

- Spondylitis
- Spondylosis
- Rheumatoid arthritis
- Osteoarthritis

Note: Made to measure service available, please provide neck circumference and front depth.

Ordering Information

Size	Neck Circ. (cm)	Front Depth (cm)	Product Code
Small	38 - 42	6.5	SFC/S
Medium	43 - 46	7	SFC/M
Large	47 - 51	8	SFC/L

Need a different size to fit your patient?

Made to Measure service is available for this product!
To order, please provide neck circumference and front depth.

