

## Push Ortho Thumb Brace CMC



### Features and Benefits

- Hand remains stable and allows good grip function
- Provides pain relief when experiencing thumb stress
- Surrounding joints remain free and doesn't limit their movement
- Pre-formed flexible plastic enclosure
- Aluminium splint can be manipulated for individual thumb contours and optimal stabilisation of the joint
- Two hook and loop fastening straps fasten around the back of the hand which allows the brace to be easily fitted by the individual
- The material and rounded edges allows for maximum comfort
- Antimicrobial
- Can be used in water and can be dried easily

### Indications for Use

- Osteoarthritis of the CMC-1 joint
- Postoperative after care of the CMC-1 joint
- Status after arthroplasty of the CMC-1 joint
- Instability of the CMC-1 joint

Measure the circumference of the affected hand and use the chart to pick the correct size. If you are in between sizes, select the larger size.



### Ordering Information

Size	MCP Circ. (cm)	Product Code	
		Left	Right
0	15 - 17.5	PUS122B	PUS121B
1	17.5 - 19.5	PUS067B	PUS064B
2	19.5 - 22.5	PUS068B	PUS065B
3	22.5 - 27	PUS069B	PUS066B

## Push Care Wrist Brace



### Features and Benefits

- Provides pressure around the wrist joint to improve stability
- Allows for adjustable pressure for individual needs, the wrist joint is not immobilised and functionality of the hand is retained
- The palm of the hand maintains freedom, ensuring a good grasp can be achieved
- Padded material around the ulnar area of the wrist
- Sympress material inner offers maximum comfort and keeps the skin dry under the brace

### Indications for Use

- Residual instability
- Moderate to serious distortions
- Strain complaints in the wrist joint
- Therapy resistant strain complaints to the wrist extensors
- Wrist fractures without dislocations
- Aftercare for wrist fractures
- Arthrosis with distension
- RSI (repetitive strain injury)



### Ordering Information

Size	Wrist Circ. (cm)	Product Code	
		Left	Right
1	13 - 15	PUS074B	PUS070B
2	15 - 17	PUS075B	PUS071B
3	17 - 19	PUS076B	PUS072B
4	19 - 21	PUS077B	PUS073B