

Push Med Wrist Brace



Features and Benefits

- Provides high levels of support and protection to the wrist joint without restriction of the functionality of the hand
- Anatomically shaped stiffeners limit palmary and dorsal flexion without inhibiting it completely
- Elastic strap provides adjustable compression
- Easy to fit
- Can be applied using one hand meaning the individual can fit it themselves
- Padded material around vulnerable areas of the wrist

Indications for Use

- Residual instability
- Moderate to serious distortions
- Strain complaints in the wrist
- Therapy-resistant strain complaints of the wrist extensors
- Wrist fractures without dislocations
- Follow-up treatment of wrist fractures
- Arthrosis with hydrops
- Repetitive strain injury



Ordering Information

Size	Wrist Circ. (cm)	Product Code	
		Left	Right
1	13 - 15	PUS082B	PUS078B
2	15 - 17	PUS083B	PUS079B
3	17 - 19	PUS084B	PUS080B
4	19 - 21	PUS085B	PUS081B

Short Wrist Brace



Features and Benefits

- Standard single compression one-piece wrist brace provides good conformity around the wrist and hand
- Removable, mouldable metal stay
- Made from Neoprene

Indications for Use

- Flaccid drop wrist
- Repetitive strain injury
- Rheumatoid arthritis
- Carpal tunnel syndrome



Ordering Information

Size	Wrist Circ. (cm)	Product Code	
		Left	Right
X-Small	10 - 12.5	SWB/E/XS/L	SWB/E/XS/R
Small	12.5 - 15	SWB/E/S/L	SWB/E/S/R
Medium	15 - 17.5	SWB/E/M/L	SWB/E/M/R
Large	17.5 - 20	SWB/E/L/L	SWB/E/L/R
X-Large	20 - 22.5	SWB/E/XL/L	SWB/E/XL/R