

Leada Multi-Use Boot



Features and Benefits

- Utilises Flex-E-Core technology for adjustment without the need for tools or heat

Indications for Use

- Plantar fasciitis
- Achilles tendonitis
- Night splinting
- Foot drop

Note: Paediatric sizes also available, see **Page 233**.

Ordering Information

Size	Max Foot Length (cm)	Brace Height (cm)	Product Code
Adult	25	35	LMU-A

Repose Wedge



Features and Benefits

- The Repose Wedge has been designed primarily for off-loading pressure from the heels and can be used wherever support or elevation of the lower legs is required. The Repose Wedge offers the same pressure reducing properties as other Repose products. Particularly suited to leg elevation to aid vascular conditions and a useful aid in posture support

Indications for Use

- The Repose Wedge is a versatile product suitable for a wide range of applications within pressure area

Ordering Information

Size	Description	Product Code
Universal	Wedge and pump	6711100
Universal	Wedge with PU strap and pump	6711111

Podus Boot



Features and Benefits

- A lightweight, quality plastic shell with poly/pile liner that helps treat and prevent lower extremity disorders associated with trauma or immobility
- An adjustable toe extension helps protect toes from bedding and the anti-rotation bar may be repositioned to assist with internal or external rotation
- Removable transfer pad on the dorsal side for limited patient ambulating
- Fits right or left foot
- Latex free

Indications for Use

- Pressure neurosis
- Ankle contractures
- Post-op hip abduction/rotation

Ordering Information

Size	Mens Shoe Size (UK)	Womens Shoe Size (UK)	Product Code
Regular	<10	<11	POD-R
X-Large	>10.5	>11.5	POD-XL

Liners (3 Pack)	
Size	Product Code
Regular	79-90560
X-Large	79-90561

More Information
Wear and Care Instructions