



PediSupport.night

Foot Postural Support Orthosis
(Night Stabilisation Splint) from
Streifeneder

The Pedi Support Night provides optimal immobilisation of the lower leg by reducing a plantar flexed position of the foot and preventing contracture of the achillies tendon.

This device maintains a plantar grade position, featuring a toe wedge to increase the plantar fascia stretch. Additionally, this product features a calf pad to offload the heel for patients in a supine position.

Indications for use:

- Plantar fasciitis
- Post-operative immobilisation of the lower leg
- Dorsal flexor - resp. peroneal paralysis
- For night postural support
- Pes equinus position

Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Product Code
Small	3.5 - 5.5	36 - 39	190R1/S
Medium	5.5 - 7	39 - 41	190R1/M
Large	7 - 9.5	41 - 44	190R1/L
X-Large	9.5 - 11	44 - 46	190R1/XL

Steeper Group, Unit 3 Stourton Link, Intermezzo Drive
Leeds LS10 1DF
Tel: +44 (0) 870 240 4133
Email: customerservices@steepergroup.com

www.steepergroup.com

STPOR171

 **steeper**
Creating life's turning points, together