

Prosthetic Limb Protection

Pro Armour Lite

PRO ARMOUR®

The Pro Armour Lite is specifically designed for patients who like to explore off the beaten track. Users can carry out their everyday activities with confidence, knowing their prosthetic foot and components have an extra layer of protection from sand, dust, dirt and splashes. Promising less hassle and more freedom, the Pro Armour Lite is simple to use, taking just seconds to don and doff, and it can save the time and hassle associated with prosthetic maintenance.

Key Features and Benefits

Splash, dust, sand, and dirt proof

The Pro Armour Lite gives users the confidence and freedom to do more, without worrying about foreign objects damaging expensive prosthetic components.

Washable and reusable

Simply hand wash with a mild detergent and air dry away from direct heat sources.

A hygienic solution

Designed to reduce the need for constant cleaning of components. This simple accessory can save both patients and prosthetists valuable time by reducing build up within the foot shell.

Lightweight construction

Made with a lightweight neoprene that weighs between 75g and 85g (depending on the size).

Multiple sizes available

The Pro Armour Lite is available in three sizes. Custom options can be requested.

Fast and simple to take on and off

The Pro Armour Lite is donned over the foot shell, using the knitted sock supplied to reduce friction where necessary. Footwear must always be worn when using the Pro Armour Lite.

Ideal for all K1-K4 patients who:

- Enjoy recreational activities such as walking, golf, fishing, cycling, camping, watersports, running, hiking and motorsports.
- Participate in domestic DIY activities including mowing lawns, gardening, washing vehicles, home maintenance.
- Work in or frequent environments where their prosthetic foot/ankle may be exposed to debris or splashes.

Foot Compatibility

The Pro Armour Lite is generally suitable for SACH, uniaxial, multiaxial, energy storing and patient adjustable feet.*

Technical Information

Pro Armour Lite					
Product Code	Size	Foot Shell Size	Weight	Material	Colour
NPALS	Small	23-24cm	75g	Neoprene	Black
NPALM	Medium	25-26cm	80g	Neoprene	Black
NPALL	Large	27-29cm	85g	Neoprene	Black
NPALC	Custom	Any shape	-	Neoprene	Any colour

*Note: The Pro Armour Lite may not fit prosthetic limbs with cosmetic covers.

