




Aluminium 30mm Structural Components

Part Number	Manufacturer	Description	Weight Limit
1* NAAA210	TL	4-Hole Pyramid Adapter	160kg
2* NAAA211	TL	4-Hole Female Adapter	160kg
3* NAAA212	TL	4-Hole Rotatable Adapter - Note: This part number is for the 4 hole base plate only. It will take either NAAA214-01 or NAAA216-01, these must be ordered separately.	160kg
4* NAAA214-01	TL	Rotatable Female Insert	160kg
5* NAAA216-01	TL	Rotatable Male Pyramid Insert - Note: This part number is for the insert only. It can be used in either the 4 bolt adapter NAAA212 or the 3-Prong laminating adapter NAAASS237, these must be ordered separately.	160kg
6* NAAA209	TL	4-Hole Single Offset Pyramid Adapter - Offset by 5mm in one direction only.	160kg
7* NAAA208	TL	4-Hole Double Offset Pyramid Adapter - Offset by 5mm in two directions.	160kg
8* NAAA100	TL	Adjustable Height Clamp Adapter - This clamp allows up to 25mm adjustment in height, with the two bolt clamp being sufficient to hold the position, unless there is likely to be any significant impact in the patient's gait, such as jumping or running.	160kg
9* NAAA135	TL	Pylon with Integral Receiver - This 30mm pylon with integral female receiver can be used with either the NAAA100 or NAAA225 tube clamps, but at 420mm long, is suitable for most transtibial applications and can be cut down to suit.	160kg

*Note: The weight limit is 123kg for the NAAA aluminium components above for patients of K4 activity level.

Manufacturer: TL = Trulife



Aluminium 30mm Structural Components

Part Number	Manufacturer	Description	Weight Limit
1* NAAA221	TL	35mm Double-Ended Adapter	160kg
2* NAAA222	TL	45mm Double-Ended Adapter	160kg
3* NAAA223	TL	60mm Double-Ended Adapter	160kg
4* NAAA224	TL	75mm Double-Ended Adapter	160kg
5* NAAA200	TL	420mm Bonded Pylon	160kg
6* NAAA202	TL	420mm Aluminium Pylon Bonded to Male Adapter	160kg
7* NAAA225	TL	30mm Clamp Adapter	160kg
8* NAAA285	TL	Foot Pyramid - This foot adapter is specifically for use with the Seattle Lightfoot, Natural or SACH, though these now come with an adapter already fitted.	136kg

Stainless Steel 30mm Structural Components

Part Number	Manufacturer	Description	Weight Limit
1* NAAASS237	TL	3-Prong Laminating Piece Stainless Steel - This part number is only for the adapter and can take either NAAA214-01 or NAAA216-01, or any component with an M36 threaded top.	136kg

*Note: The weight limit is 123kg for the NAAA aluminium components above for patients of K4 activity level.

Manufacturer: TL = Trulife