



Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

Stainless Steel College Park Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
1	N4HPA-S	CP	4-Hole Pyramid Adapter, Standard	181kg
2	N4HPR-S	CP	4-Hole Pyramid Receiver, Standard	181kg
3	NTA-S250-30	CP	Low-Profile Tube Adapter, 250mm length, 30mm diameter	136kg
4	NTA-S430-30	CP	Low-Profile Tube Adapter, 430mm length, 30mm diameter	136kg
5	NTC-S-30	CP	Tube Clamp, 30mm diameter	136kg
6	NSS-TTK-PA	CP	Transtibial Kit with Pyramid Adapter (30mm Tube Adapter - 250mm length, 30mm Tube Clamp, 4-Hole Pyramid Adapter)	136kg
7	NSS-TFK-PA	CP	Transfemoral Kit with Pyramid Adapter (30mm Tube Adapter - 430mm length, 30mm Tube Clamp, 4-Hole Pyramid Adapter)	136kg

Manufacturer: CP = College Park



Weight Limit Note: Please note that when ordering components from different ranges for the same limb build, the lowest rated component becomes the maximum weight limit.

34mm Heavy Duty Stainless Steel College Park Structural Components				
Part Number	Manufacturer	Description	Weight Limit	
1	N4HPA-S	CP	4-Hole Pyramid Adapter, Standard	181kg
2	N4HPR-S	CP	4-Hole Pyramid Receiver, Standard	181kg
3	NTA-S250-34	CP	Low-Profile Tube Adapter, 250mm length, 34mm diameter	181kg
4	NTA-S430-34	CP	Low-Profile Tube Adapter, 430mm length, 34mm diameter	181kg
5	NTC-S-34	CP	Tube Clamp, 34mm diameter	181kg
6	NHD-TTK-PA	CP	Transtibial Heavy Duty Kit with Pyramid Adapter (34mm Tube Adapter - 250mm length, 34mm Tube Clamp, 4-Hole Pyramid Adapter)	181kg
7	NHD-TFK-PA	CP	Transfemoral Heavy Duty Kit with Pyramid Adapter (34mm Tube Adapter - 430mm length, 34mm Tube Clamp, 4-Hole Pyramid Adapter)	181kg

Manufacturer: CP = College Park