

# ContexGel.sleeve

The Streifeneder ContexGel.sleeves ensure correct suspension and fit in lower and upper limb prosthetic applications. They're effective at keeping air out of elevated system sockets and keeping the prosthesis firmly in place. Smaller sizes are suitable for upper limb below-elbow prosthetic sockets; larger sizes are suitable for below-knee amputees.

The latest bi-elastic knit design prevents tension and pinching on the knee joint, with 3mm antibacterial Polymer Gel for optimal residual limb health.

## Features and Benefits

- Conical shaped sleeve
- Available in two colours: mixed anthracite or sand
- Round-knit design with no lateral seams
- Soft, abrasion-resistant textile cover which can be easily donned or doffed
- Features an outer matrix to eliminate elongation in the distal and proximal area
- Excellent circumferential elasticity compensates for residual limb volume changes



# Liners - Streifeneder - ContexGel.sleeve

## Ordering Information

| Colour: Mixed Anthracite |           |           |        |               |
|--------------------------|-----------|-----------|--------|---------------|
| Size                     | Distal    | Proximal  | Length | Part Number   |
| S                        | 15 - 23cm | 23 - 35cm | 37cm   | NMA30L10/S    |
| M                        | 21 - 31cm | 29 - 44cm | 37cm   | NMA30L10/M    |
| L                        | 28 - 43cm | 36 - 52cm | 38.5cm | NMA30L10/L    |
| XL                       | 34 - 48cm | 42 - 58cm | 38.5cm | NMA30L10/XL   |
| XXL                      | 38 - 52cm | 43 - 64cm | 38.5cm | NMA30L10/XXL  |
| XXXL                     | 42 - 56cm | 50 - 70cm | 38.5cm | NMA30L10/XXXL |

| Colour: Sand |           |           |        |               |
|--------------|-----------|-----------|--------|---------------|
| Size         | Distal    | Proximal  | Length | Part Number   |
| S            | 15 - 23cm | 23 - 35cm | 37cm   | NMA30L50/S    |
| M            | 21 - 31cm | 29 - 44cm | 37cm   | NMA30L50/M    |
| L            | 28 - 43cm | 36 - 52cm | 38.5cm | NMA30L50/L    |
| XL           | 34 - 48cm | 42 - 58cm | 38.5cm | NMA30L50/XL   |
| XXL          | 38 - 52cm | 43 - 64cm | 38.5cm | NMA30L50/XXL  |
| XXXL         | 42 - 56cm | 50 - 70cm | 38.5cm | NMA30L10/XXXL |

## Size Information

ContexGel Knee Sleeves are available in six sizes. When determining patient size, please position the sleeve so that the middle section is sat over the knee/elbow joint.

To determine the optimal size of the knee sleeve, measure the distal socket circumference and the proximal thigh circumference (approx. 15cm above the centre of the knee joint).

## ContexGel.protect



To avoid premature tearing of the ContexGel knee sleeves, we recommend using the ContexGel.protect as a protective cover for sharp-edged prosthetic sockets. This protective cover is made of durable and expandable fabric and has an inner silicone edge on the distal side to prevent slipping.

| ContexGel.protect cover |           |                 |        |             |
|-------------------------|-----------|-----------------|--------|-------------|
| Size                    | Side      | For Knee Sleeve | Length | Part Number |
| S                       | Universal | Size S          | 19cm   | NMA3S50/S   |
| M                       | Universal | Size M          | 19cm   | NMA3S50/M   |
| L                       | Universal | Size L          | 22cm   | NMA3S50/L   |
| XL                      | Universal | Size XL & XXL   | 22cm   | NMA3S50/XL  |

## Protective Inner Sleeve



These sleeve protectors from Össur are made of durable fabric that extends the life of the suspension sleeve and protects it from damage caused by the socket's edge. Silicone pimples on the bottom edge of the sleeve hold it securely in place on the socket.

| Protective Sleeve |             |
|-------------------|-------------|
| Size              | Part Number |
| Small             | NM000202    |
| Medium            | NM000203    |
| Large             | NM000204    |

## Technical Information

|                  |             |
|------------------|-------------|
| Material         | Polymer Gel |
| Length           | 37 - 38.5cm |
| Amputation Level | Transtibial |
| Warranty         | 3 months    |