

Nylon Protective Sleeves

Fine, highly elastic nylon protective sleeves that act as a second skin under the cotton residual limb sock.

The material textures ensure adherence of the sheath to the skin, and guarantees that the cotton residual limb sock dons easily.

Should be cared for with lukewarm water and a mild neutral detergent.

Features and Benefits

- High level of wearing comfort
- Colour: white
- 2 sheaths in each pack
- Washable



Socks - Uniprox - Nylon Protective Sleeves

Ordering Information

Part Number	Proximal Width	Distal Width	Length	Colour
NG-910-B-30	16cm	8cm	30cm	White
NG-910-B-35	16cm	8cm	35cm	White
NG-910-B-40	16cm	8cm	40cm	White
NG-910-B-45	16cm	8cm	45cm	White
NG-910-B-50	16cm	8cm	50cm	White
NG-910-B-55	16cm	8cm	55cm	White
NG-910-B-60	16cm	8cm	60cm	White
NG-910-C-30	18cm	10cm	30cm	White
NG-910-C-35	18cm	10cm	35cm	White
NG-910-C-40	18cm	10cm	40cm	White
NG-910-C-45	18cm	10cm	45cm	White
NG-910-C-50	18cm	10cm	50cm	White
NG-910-C-55	18cm	10cm	55cm	White
NG-910-C-60	18cm	10cm	60cm	White
NG-910-D-30	20cm	10cm	30cm	White
NG-910-D-35	20cm	10cm	35cm	White
NG-910-D-40	20cm	10cm	40cm	White
NG-910-D-45	20cm	10cm	45cm	White
NG-910-D-50	20cm	10cm	50cm	White
NG-910-D-55	20cm	10cm	55cm	White
NG-910-D-60	20cm	10cm	60cm	White