



Durable, yet lightweight.

Designed with the basic principles of durability and ease of use firmly in mind.

Using high strength aerospace alloy inside a flexible PVC shell, these prosthetic hands are also compatible with our range of full-length Elegance™ cosmetic gloves in 19 lifelike TrueFinish™ shades for enhanced realism.

Our advanced materials and manufacturing processes ensure that these cable and spring operated prosthetic hands are detailed, durable and lightweight options for patients.

Featuring:

- Body-powered cable control or manual operation
- Voluntary opening - front pull, back pull or spring close options
- 3 cable terminations supplied
- Four sizes available from child to adult
- Plain hand base, M12 x 1.5 or ½ x 20 TPI threaded stem



Cable and Spring Operated Mechanical Hands



To Order

Prefix the number shown in the table with TG for TrueFinish PVC or TSG for TrueFinish Silicone.

Then, add the preferred skin shade. E.g. **TSG300/E4** orders a 8¼ left adult silicone glove in TrueFinish colour E4.

To order a colour swatch use part number **E24655**.



| Ordering Information - Spring Operated Hands | | | | | | | | |
|----------------------------------------------|-------------------------|--------------------------|------|--------------|---------------------|----------------|--------------------------------------------------|--------------|
| Part No. Plain Base | Part No. M12 x 1.5 Stud | Part No. ½ x 20 TPI Stud | Size | Left / Right | Application | Wrist Diameter | TG = TrueFinish PVC TSG = TrueFinish Silicone | |
| | | | | | | | Male Glove | Female Glove |
| SPA675L | SPA675LEU | SPA675LNA | 6 ¾ | Left | Child (Approx.8-12) | 45mm | 310 | 310 |
| SPA675R | SPA675REU | SPA675RNA | 6 ¾ | Right | Child (Approx.8-12) | 45mm | 311 | 311 |
| SPA725L | SPA725LEU | SPA725LNA | 7 ¼ | Left | Adolescent / Adult | 50mm | 306 | 304 |
| SPA725R | SPA725REU | SPA725RNA | 7 ¼ | Right | Adolescent / Adult | 50mm | 307 | 305 |
| SPA775L | SPA775LEU | SPA775LNA | 7 ¾ | Left | Adult | 54mm | 302 | 312 |
| SPA775R | SPA775REU | SPA775RNA | 7 ¾ | Right | Adult | 54mm | 303 | 313 |
| SPA825L | SPA825LEU | SPA825LNA | 8 ¼ | Left | Adult | 54mm | 300 | |
| SPA825R | SPA825REU | SPA825RNA | 8 ¼ | Right | Adult | 54mm | 301 | |

| Ordering Information - Cable Operated Hands - Front Pull Voluntary Opening | | | | | | | | |
|----------------------------------------------------------------------------|-------------------------|--------------------------|------|--------------|---------------------|----------------|--------------------------------------------------|--------------|
| Part No. Plain Base | Part No. M12 x 1.5 Stud | Part No. ½ x 20 TPI Stud | Size | Left / Right | Application | Wrist Diameter | TG = TrueFinish PVC TSG = TrueFinish Silicone | |
| | | | | | | | Male Glove | Female Glove |
| MHA675LF | MHA675LFEU | MHA675LFNA | 6 ¾ | Left | Child (Approx.8-12) | 45mm | 310 | 310 |
| MHA675RF | MHA675RFEU | MHA675RFNA | 6 ¾ | Right | Child (Approx.8-12) | 45mm | 311 | 311 |
| MHA725LF | MHA725LFEU | MHA725LFNA | 7 ¼ | Left | Adolescent / Adult | 50mm | 306 | 304 |
| MHA725RF | MHA725RFEU | MHA725RFNA | 7 ¼ | Right | Adolescent / Adult | 50mm | 307 | 305 |
| MHA775LF | MHA775LFEU | MHA775LFNA | 7 ¾ | Left | Adult | 54mm | 302 | 312 |
| MHA775RF | MHA775RFEU | MHA775RFNA | 7 ¾ | Right | Adult | 54mm | 303 | 313 |
| MHA825LF | MHA825LFEU | MHA825LFNA | 8 ¼ | Left | Adult | 54mm | 300 | |
| MHA825RF | MHA825RFEU | MHA825RFNA | 8 ¼ | Right | Adult | 54mm | 301 | |

| Ordering Information - Cable Operated Hands - Back Pull Voluntary Opening | | | | | | | | |
|---------------------------------------------------------------------------|-------------------------|--------------------------|------|--------------|---------------------|----------------|--------------------------------------------------|--------------|
| Part No. Plain Base | Part No. M12 x 1.5 Stud | Part No. ½ x 20 TPI Stud | Size | Left / Right | Application | Wrist Diameter | TG = TrueFinish PVC TSG = TrueFinish Silicone | |
| | | | | | | | Male Glove | Female Glove |
| MHA675LB | MHA675LBEU | MHA675LBNA | 6 ¾ | Left | Child (Approx.8-12) | 45mm | 310 | 310 |
| MHA675RB | MHA675RBEU | MHA675RBNA | 6 ¾ | Right | Child (Approx.8-12) | 45mm | 311 | 311 |
| MHA725LB | MHA725LBEU | MHA725LBNA | 7 ¼ | Left | Adolescent / Adult | 50mm | 306 | 304 |
| MHA725RB | MHA725RBEU | MHA725RBNA | 7 ¼ | Right | Adolescent / Adult | 50mm | 307 | 305 |
| MHA775LB | MHA775LBEU | MHA775LBNA | 7 ¾ | Left | Adult | 54mm | 302 | 312 |
| MHA775RB | MHA775RBEU | MHA775RBNA | 7 ¾ | Right | Adult | 54mm | 303 | 313 |
| MHA825LB | MHA825LBEU | MHA825LBNA | 8 ¼ | Left | Adult | 54mm | 300 | |
| MHA825RB | MHA825RBEU | MHA825RBNA | 8 ¼ | Right | Adult | 54mm | 301 | |

Steeper Group, Unit 3 Stourton Link, Intermezzo Drive, Leeds, LS10 1DF

www.steepergroup.com

STPPR038 Issue 1 2017