



Walking Stick Holder

The walking stick holder is connected to the prosthesis with an angled tube. Supplied with 2 1/8" adapter plate. Walking stick not included.

Part Number	Description	Left/Right	Thread	Length	Weight
B05640	Walking Stick Holder - Adult	Left or Right	3/8" BSW	180mm (7 1/8")	146g (5 oz)



Plain Hook

For use in manual handling and lifting applications.

Part Number	Description	Left/Right	Thread	Length	Weight
B05138	Plain Hook - Adult Short	Left or Right	3/8" BSW	148mm (6")	109g (3 4/5 oz)
B05140	Plain Hook - Adult Long	Left or Right	3/8" BSW	312mm (12 1/4")	198g (7 oz)



'C' Hook

A multi-purpose hook for use with gardening and farming tools (e.g. broom, shovel, fork, wheelbarrow) where pushing and pulling are required. The head can swivel and flex, providing a more adaptive device. The plain 'C' Hook retains functional head shape but is attached to the stem in a fixed position.

Part Number	Description	Left/Right	Thread	Length	Weight
B00551	'C' Hook - Adult	Left or Right	3/8" BSW	141mm (5 1/2")	190g (6 1/2 oz)
B00770	Plain 'C' Hook - Adult	Left or Right	3/8" BSW	152mm (6")	132g (4 1/2 oz)

Note: A full range of adapter plates are available. Please see page 138. Thread locking compound must be used when assembling adapters to functional devices.

Functional Devices - Passive Functional Devices



Potato Holder

Allows items such as vegetables, fish, bread and meat etc, to be held securely during food preparation, while leaving the sound hand free to peel or cut.

Part Number	Description	Left/Right	Thread	Length	Weight
B00538	Potato Holder - Adult	Left or Right	$\frac{3}{8}$ " BSW	130mm (5 1/8")	100g (3 1/2 oz)



Vegetable Peeler

Sharp cutting edge allows for easy peeling of fruit or vegetables.

Part Number	Description	Left/Right	Thread	Length	Weight
B11476	Vegetable Peeler - Adult	Left or Right	$\frac{3}{8}$ " BSW	125mm (5")	16g (1/2 oz)

Note: A full range of adapter plates are available. Please see page 138.

Thread locking compound must be used when assembling adapters to functional devices.



Table Fork

Large standard fork manufactured using stainless steel.

Part Number	Description	Left/Right	Thread	Length	Weight
1	B05355	Table Fork - Adult	Left or Right	3/8" BSW	166mm (6 1/2") 42g (1 1/2 oz)
1a	NGO3F83L	Angled Table Fork - Adult	Left	3/8" BSW	170mm (6 3/4") 64g (2 1/4 oz)
	NGO3F83R	Angled Table Fork - Adult	Right	3/8" BSW	170mm (6 3/4") 64g (2 1/4 oz)



Table Knife

Large standard knife manufactured using stainless steel.

Part Number	Description	Left/Right	Thread	Length	Weight
2	B05354	Table Knife - Adult	Left or Right	3/8" BSW	173mm (7") 48g (1 1/2 oz)
2a	NGO3F81	Steak Knife - Adult	Left or Right	3/8" BSW	130mm (5 1/8") 40g (1 1/2 oz)



Table Spoon

Large standard spoon manufactured using stainless steel.

Part Number	Description	Left/Right	Thread	Length	Weight
3	B05356	Table Spoon - Adult	Left or Right	3/8" BSW	140mm (5 1/2") 42g (1 1/2 oz)
3a	NGO3F82L	Angled Table Spoon - Adult	Left	3/8" BSW	175mm (6 7/8") 82g (3 oz)
	NGO3F82R	Angled Table Spoon - Adult	Right	3/8" BSW	175mm (6 7/8") 82g (3 oz)

Note: A full range of adapter plates are available. Please see page 138. Thread locking compound must be used when assembling adapters to functional devices.

Functional Devices - Passive Functional Devices



Fishing Rod Holder

Used to support a fishing rod, allowing the sound hand freedom to control the reel.

Part Number	Description	Left/Right	Thread	Length	Weight
B05882	Fishing Rod Holder	Left or Right	3/8" BSW	155mm (6")	143g (5 oz)



Typing Appliance

Supplied with an angled head and rubber tip for operation of a computer keyboard.

Part Number	Description	Left/Right	Thread	Length	Weight
B05362	Typing Appliance - Adult	Left or Right	3/8" BSW	115mm (4 1/2")	48g (1 1/2 oz)
B05481	Replacement Rubber Tip	N/A	N/A	N/A	N/A



Cue Rest

Designed to support the cue, this device has a rotating head to provide a stable base and three different cueing heights.

Part Number	Description	Left/Right	Thread	Length	Weight
B07318	Cue Rest - Adult	Left or Right	3/8" BSW	174mm (7")	56g (2 oz)

Note: A full range of adapter plates are available. Please see page 138.
Thread locking compound must be used when assembling adapters to functional devices.



Thread Conversion Adapter

Can be applied to convert external UK thread to US or EU thread, available in three sizes.

Part Number	Description	Supplied
B06005	Conversion Sleeve - External Thread 1/2" UNF x 20 TPI, Internal Thread 3/8" BSW To convert UK External Thread to US External Thread	Each
B21871	Conversion Sleeve - External Thread M12 x 1.5mm, Internal Thread 3/8" BSW To convert UK External Thread to European External Thread	Each
B23372	Conversion Sleeve - External Thread 1/2" UNF x 20 TPI, Internal Thread 5/16" BSF To convert Child UK External Thread to US External Thread	Each



Extension Bar

An extension piece which may be used to increase length between two prosthetic components.

Part Number	Description	Left/Right	Thread	Weight
B00445	Extension Bar 3" - Adult	N/A	3/8" BSW External and Internal	71g (2 1/2 oz)
B00446	Extension Bar 4" - Adult	N/A	3/8" BSW External and Internal	96g (3 1/2 oz)

Note: A full range of adapter plates are available. Please see page 138.
Thread locking compound must be used when assembling adapters to functional devices.

Functional Devices - Passive Functional Devices

Upper Limb
Prosthetics

Ordering
Information

Myoelectric
Hands

Myoelectric
Elbows

Myoelectric
Power & Inputs

Cable Operated
& Spring Hands

Functional
Devices

Cosmetic Hands
& Gloves

Modular Arm
Systems

Wrist Units

Mech. Elbow
Units & Joints

Suspension &
Control Systems


Materials &
Tools



Nail Brush Holder

This stainless steel adapter has two drilled holes for attachment to the object of a user's choice, e.g. nail brush, wire brush, files, paint roller etc. Supplied as shown.

Part Number	Description	Left/Right	Thread	Length	Weight
B05353	Nail Brush Holder - Adult	Left or Right	3/8" BSW	114mm (4 1/2")	54g (2 oz)



Telescopic Dressing Stick

Extendible device designed to facilitate buttoning and unzipping using either mouth or residual digits. May be used with a mouthpiece.

Part Number	Description	Left/Right	Length	Weight
1	B05553 Telescopic Dressing Appliance	N/A	Max 340mm / Min 220mm	65g (2 1/2 oz)
1a	SB12202 Mouthpiece	N/A		15g (1/2 oz)

Note: A full range of adapter plates are available. Please see page 138.

Thread locking compound must be used when assembling adapters to functional devices.



Motoring Ball and Clamp - Adult*

For driving applications. A rim-mounted single clamp arrangement fits most types of steering wheel.

***Available in the UK only.**

Note: E22252A is supplied with link bar to increase stability and can be used with vehicles fitted with air bags.

Part Number	Description	Left/Right	Weight
E22252	Motoring Ball and Clamp - Adult	N/A	349g (12 1/3 oz)
E22252A	Motoring Clamp with Link - Adult	N/A	255g (9 oz)



Ball and Clamp

This device may be used as a handlebar-mounted aid for adults and children wishing to cycle.

Note: Not for use with motorised vehicles.

Part Number	Description	Left/Right	Weight
SP118	Ball and Clamp	N/A	115g (4 oz)

Note: A full range of adapter plates are available. Please see page 138. Thread locking compound must be used when assembling adapters to functional devices.

Driving Cup*

For driving applications, this device is screwed into a palm mounted fitting (B05718). It must be used with an appropriate steering ball and clamp. A fitted PVC cover gives a secure level of attachment while allowing rapid removal if necessary.



Note: Alternative cup sizes may be manufactured. Please contact Customer Services with your requirements.
*Available in the UK only.

Part Number	Description	Left/Right	Thread	Weight
1 B00447	Driving Cup - Adult	N/A	3/8" BSW	126g (4 1/2oz)
1a B05718	Palm Mounting Fitting	N/A	N/A	26.4g (7/8 oz)



Driving Cup on Stem*

For driving applications. Must be used with an appropriate steering ball and clamp. A fitted PVC cover gives a secure level of attachment while allowing rapid removal if necessary.

Note: Alternative cup sizes may be manufactured. Please contact Customer Services with your requirements.
*Available in the UK only.

Part Number	Description	Left/Right	Thread	Weight
B00451	Driving Cup on Stem - Adult	N/A	3/8" BSW	133g (4 1/2oz)

Note: A full range of adapter plates are available. Please see page 138.
Thread locking compound must be used when assembling adapters to functional devices.