

COLLEGE PARK - ORDER FORM

Fax to RSLSTEEPER Sales Desk: Fax 0870 240 4779 Tel 0870 240 4133

RSLSTEEPER, Leeds Manufacturing Centre, Unit 7, Hunslet Trading Estate, Severn Road, Leeds LS10 1BL

PATIENT INFORMATION

PATIENT'S NAME:

SEX: M F CAUCASIAN BROWN





AGE: WEIGHT: HEIGHT:

AMPUTATION LEVEL (TF, TT, KD, HD...): LEFT RIGHT SIZE (Cm):

PROSTHETIST: BRANCH / CONTRACTOR:

FOOT STYLE

Patient Impact level: Please circle **1** **2** **3** (See below for Descriptions of Impact Levels)

| | | | |
|---|---|---|---|
|  <p>TRÉS <input type="checkbox"/></p> <p>Multi Axial/ Dynamic/ non-adjustable foot for low to moderate activity. Up to level 2 only</p> |  <p>TRIBUTE <input type="checkbox"/></p> <p>Multi Axial/ Dynamic/ adjustable foot for moderate activity. Up to level 2 only</p> |  <p>VENTURE <input type="checkbox"/></p> <p>Multi Axial/ Dynamic foot for any activity level.</p> |  <p>TRUSTEP <input type="checkbox"/></p> <p>Polycentric/ Multi Axial/ Dynamic foot for any activity level.</p> |
|---|---|---|---|

ATTACHMENTS

| | | | |
|--|---|---|---|
| <p>Pyramid <input type="checkbox"/></p> <p>125kg (sizes 25-30) weight limit 100kg (sizes 21-30) weight limit</p> | <p>Pyramid <input type="checkbox"/></p> <p>100kg weight limit</p> | <p>Pyramid <input type="checkbox"/></p> <p>Alignable Exo <input type="checkbox"/></p> <p>125kg weight limit</p> | <p>Pyramid <input type="checkbox"/></p> <p>136kg weight limit</p> <p>Alignable Exo <input type="checkbox"/></p> <p>160kg weight limit</p> |
|--|---|---|---|

PAEDIATRIC FOOT

TRUPER  Pyramid Alignable Exo Growth Plate Kit

Impact Level Descriptions

- 1. Low** : Daily activities include mostly level ground walking, moving around in the home and the community
- 2. Moderate** : Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (This does not include running)
- 3. High** : Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports

| Part Number | Description | Quantity |
|-------------|-------------|----------|
| | | |