

# THE DEFENDER PAEDIATRIC FOOT

TAKING A STAND FOR CHILDREN EVERYWHERE.

FREEDOM TO RUN.

FREEDOM TO PLAY.

FREEDOM TO JUMP.

FREEDOM TO LIVE.



# FREEDOM SERIES™

Finally. Performance and Value  
in Paediatric Prosthetics.

Freedom Innovations, the leader in advanced technology prosthetic feet, now offers children the freedom and flexibility to live life to its fullest.

The **Defender** Paediatric foot is newly designed to ensure smooth rollover and exceptional ground compliance. The best way to effectively increase energy storage and return is to increase the length of the carbon fibre keel – a prominent feature of the **Defender**. And, the **Defender** offers greater energy storage and return than any paediatric foot on the market.

The **Defender's** integrated paediatric pyramid attachment and patented technology eliminates bolt holes through all areas of the foot that must flex to ensure optimal strength and performance. The foot is virtually maintenance-free – no bumpers to replace. Most importantly, Freedom's paediatric foot has been uniquely designed to optimise the varied gait of young children, yet is robust enough to endure the physical demands of a child.

## The Defender:

- Sizes 16-21 cm, six stiffness categories
- Smooth rollover
- Exceptional ground compliance
- Increased toe lever length
- Exceptional energy storage, return
- Maintenance free



## Stiffness Chart:

FOOT STIFFNESS CATEGORY	WEIGHT #
1 Paediatric	20-30
2 Paediatric	31-40
3 Paediatric	41-50
4 Paediatric	51-60
5 Paediatric	61-70
6 Paediatric	71-80
1 Adult	81-100

## ORDERING EXAMPLE - NFSP3418L

	Defender code	Category	Size	Left/Right
NFSP	3	4	18	L



Products Division, Unit 7, Hunslet Trading Estate, Severn Road, Leeds LS10 1BL, United Kingdom

Tel: +44 (0)870 240 4133 Fax: +44 (0)870 240 4779

Email: [sales@rslsteeper.co.uk](mailto:sales@rslsteeper.co.uk) Web: [www.rslsteeper.com](http://www.rslsteeper.com)

ONCE AGAIN, THE LEADER IN INNOVATION, VALUE AND TECHNOLOGY.