

# ELEGANCE

Extra Length Silicone & PVC  
Cosmetic Gloves

**AVAILABLE  
IN SEVEN  
STYLES**



## ELEGANCE Gloves

- Extra length as standard in 7 styles
- Recommended for RSL Steeper *Select*, Otto Bock™, Motion Control™ and Hosmer™ adult and youths, myoelectric and mechanical hand ranges
- A choice between the Standard colour range or the advanced micro-pigmentation Truefinish™ (available for added realism and depth)
- Both ranges available in 19 colour shades



Available in 19 colour shades, 5 sizes & 7 styles

## ELEGANCE Silicone

By using Tufsil 3™ our innovative polymer, ELEGANCE Silicone gloves offer the toughness associated with PVC but provide the ease of cleaning provided by silicones.

In addition, Tufsil 3™ allows additional stretch, reducing the energy expended by the hand motor each time it opens. Properties of this polymer include improved tear resistance proven by clinical evaluations and up to 15% stretch, to allow even bulbous sockets to be accommodated.

All ELEGANCE Silicone gloves are easy to fit and remove, allowing quick access to the hardware for maintenance and adjustment. The Silicone gloves are simple to keep clean using soap and water.

## ELEGANCE PVC

ELEGANCE PVC gloves provide a cost-effective solution for electric and mechanical hands when the extra toughness of PVC is required.

**ELEGANCE Silicone and ELEGANCE PVC** is available in the following styles:

Size	Approx Length	Application	Left Silicone	Right Silicone	Left PVC	Right PVC	Colour Ranges
<b>STANDARD COLOURS</b>							
8¼	645mm	Male	SG300	SG301	CG300	CG301	E0 to E18
7¾	645mm	Male	SG302	SG303	CG302	CG303	E0 to E18
7¾	620mm	Female	SG312	SG313	CG312	CG313	E0 to E18
7¼	645mm	Female	SG304	SG305	CG304	CG305	E0 to E18
7¼	645mm	Youths	SG306	SG307	CG306	CG307	E0 to E18
7	520mm	Female/Youths	SG308	SG309	CG308	CG309	E0 to E18
6¾	460mm	Youths	SG310	SG311	CG310	CG311	E0 to E18
<b>TRUEFINISH™ COLOURS</b>							
8¼	645mm	Male	TSG300	TSG301	TG300	TG301	E0 to E18
7¾	645mm	Male	TSG302	TSG303	TG302	TG303	E0 to E18
7¾	620mm	Female	TSG312	TSG313	TG312	TG313	E0 to E18
7¼	645mm	Female	TSG304	TSG305	TG304	TG305	E0 to E18
7¼	645mm	Youths	TSG306	TSG307	TG306	TG307	E0 to E18
7	520mm	Female/Youths	TSG308	TSG309	TG308	TG309	E0 to E18
6¾	460mm	Youths	TSG310	TSG311	TG310	TG311	E0 to E18

### ORDERING EXAMPLES

#### TSG300 / E5

Orders a 8¼ Left Silicone glove in Truefinish colour E5

#### CG300 / E5

Orders a 8¼ Left PVC glove in standard finish colour E5

#### Swatches are available on request:

E24654 E0 to E18 Standard colours

E24655 E0 to E18 Truefinish colours



#### RSLSTEEPER

Products Division

Unit 7, Hunslet Trading Estate  
Severn Road, Leeds LS10 1BL

Tel: +44(0) 870 240 4133

Fax: +44(0) 870 240 4779

Email: sales@rslsteeper.co.uk

Web: www.rslsteeper.com