

RSLSTEEPER Announces

The Highlander Max® Prosthetic Foot & Kit

► **Rated for patients weighing up to 500 lbs. (226kg)**

The new Highlander Max® Prosthetic Foot Kit provides a comprehensive solution to meet the needs of amputees weighing between 366lbs (166kg) and 500 lbs (226kg)

- Additional carbon fiber layers provide durability and energy return.
- An integral 34mm male pyramid mates with the 34mm Heavy Duty Component Kit for a simple, yet dependable system set-up.

How to Order:

- Highlander Max Prosthetic Foot
- Part No NFS3 # ***
- (# = Keel-Select from table below)
- (***) = Foot Size & Side EG 26R)
- 34mm Heavy Duty Component Kit
- Part No NAC 202 (Kit includes below)
 - Two 34mm Tube Clamps
 - One 4 Ear Connector
 - One Pylon

Keel Selection Chart				
Weight	Impact level			Weight
(lbs)	Low	Mod	High	(kg)
286-325			10	131-147
326-365		10	11	148-166
366-420	10	11	12	166-190
421-480	11	12	13	191-217
481-500	12	13		218-226

For more information, please contact RSLSTEEPER Customer Sales Desk 0870 240 4133



Warranty & Returns

The Highlander Max Prosthetic Foot is sold with a 36 month warranty. It is a custom made, non-returnable product.

The Heavy Duty Component Kit is sold with a 24 month warranty and may be returned within 60 days if unused.

Sales Desk
Unit 7, Hunslet Trading Estate
Severn Road, Leeds LS10 1BL
Phone: 0870 240 4133
Fax: 0870 240 4779
www.rslsteeper.com



