

motionSupport[®] Insoles from RSLSteeper


An innovative range of ready-made insoles featuring unique design and construction. Keeping patients active for longer.



motionSupport® Insoles from



- Full length polyurethane insoles with 3/4 length fibreglass reinforced shell offering **flexibility, durability** and **excellent value for money**
- Fibreglass extension on the lateral aspect **encourages the transition** from low to high gear propulsion in late stance phase, leading to a more **propulsive gait**
- **Slimmer profile** for **easier fitting** in footwear
- Three arch heights available in each size to **optimise sagittal plane control**
- **Flexible construction** under the first metatarsal head allows first ray to plantarflex
- A small metatarsal button works to **reduce pressure** on the metatarsal heads
- Full length construction offers **improved forefoot control**
- **Durable** alcantara top cover



Priced at only £21.99 per pair excluding VAT

Feet come in all shapes and sizes. *motionSupport*® Insoles offer increased choice and improved shaping to enhance fit and optimise control, ensuring a better outcome for the patient.

| Insole Description | Product Code |
|-----------------------------------|----------------|
| Normal Foot - Grey, size 36/37 | KRI-7052=36/37 |
| Normal Foot - Grey, size 38/39 | KRI-7052=38/39 |
| Normal Foot - Grey, size 40/41 | KRI-7052=40/41 |
| Normal Foot - Grey, size 42/43 | KRI-7052=42/43 |
| Normal Foot - Grey, size 44/45 | KRI-7052=44/45 |
| Normal Foot - Grey, size 46/47 | KRI-7052=46/47 |
| High Arch Foot - Blue, size 36/37 | KRI-7053=36/37 |
| High Arch Foot - Blue, size 38/39 | KRI-7053=38/39 |
| High Arch Foot - Blue, size 40/41 | KRI-7053=40/41 |
| High Arch Foot - Blue, size 42/43 | KRI-7053=42/43 |
| High Arch Foot - Blue, size 44/45 | KRI-7053=44/45 |
| High Arch Foot - Blue, size 46/47 | KRI-7053=46/47 |
| Low Arch Foot - Brown, size 36/37 | KRI-7055=36/37 |
| Low Arch Foot - Brown, size 38/39 | KRI-7055=38/39 |
| Low Arch Foot - Brown, size 40/41 | KRI-7055=40/41 |
| Low Arch Foot - Brown, size 42/43 | KRI-7055=42/43 |
| Low Arch Foot - Brown, size 44/45 | KRI-7055=44/45 |
| Low Arch Foot - Brown, size 46/47 | KRI-7055=46/47 |

How to order:

To order *motionSupport*® Insoles, please use the product codes on the left.

For more information or to place an order please contact RSLSteeper:

Tel: 0113 272 5688

Fax: 0113 272 5689

Email: enquiries@rslsteeper.com

