



medi<sup>®</sup>  
New product Information


medi<sup>®</sup> catalogue supplement

# Stump Socks


## medi<sup>®</sup> RELAX hip disartic/single-legged pantyhose



**NEW!**



Allivates phantom pain, sensations and idiopathic residual limb pain without side effects!



### INDICATIONS

- Protection of the residual limb
- Phantom pain and phantom sensations, idiopathic residual limb pain
- Hip disarticulation, hemipelvectomy, pseudo-hip disarticulation, short above-knee amputations

### PROPERTIES

- Relieves or eliminates phantom pain, phantom sensations and idiopathic residual limb pain due to isolation from electromagnetic effects thanks to the metallurgical knit (*Umbrellan<sup>®</sup>*)
- Knit with temperature-regulating thread (*Clima-Comfort<sup>®</sup>*):
  - warming and cooling function, depending on environmental temperature
  - highly elastic knit with SOFT topbands for comfortable, reliable fit – without slipping or constriction
  - pleasantly comfortable due to the soft, skin-friendly Umbrellan<sup>®</sup> and Clima-Comfort<sup>®</sup> material
  - no heat build-up, as stocking is breathable and moisture-regulating thanks to the well-known Clima-Comfort<sup>®</sup> technology
- Generally suitable for prostheses, but this must be checked by the practitioner in each case.



**Umbrellan<sup>®</sup> is a unique, patented knitted fabric that combines maximum elasticity with an excellent electromagnetic screen.**

The materials used on the metallurgical base are hypoallergenic, climate-regulating and have antibacterial properties.

Its functionality, therapeutic benefit and quality have been statistically significantly demonstrated.

**Prospective randomised study available as Art. No. 70.914.**

### TECHNISCHE DATEN

Material	Umbrellan <sup>®</sup> with Clima-Comfort <sup>®</sup>
Sizes / Length	see product table
Amputation level	Hip disarticulation, hemipelvectomy, pseudo-hip disarticulation, short above-knee amputations
Cleaning	Hand wash at max. 30°C using medi <sup>®</sup> care neutral soap M. 870.902
Supplied	Umbrellan <sup>®</sup> User Set (medi <sup>®</sup> pain diary, medi <sup>®</sup> pain scale)



Schadstoffgeprüfte Textilien  
Tested for harmful substances  
nach / according to  
Öko-Tex Standard 100



142a

142b



medipro®  
New product Information

medipro® catalogue supplement

# Stump Socks

## medipro® RELAX hip disartic/single-legged pantyhose



### Ordering tip:

- Use the example here on the left for your fax order!
- Delivery time approx. 10 working days


### Product information:

Each medipro® RELAX hip disartic-/single-legged pantyhose is made to measure individually for each user.

Please give all the order data to ensure an optimum fit.

### SIZE SELECTION AND INFORMATION REQUIRED FOR THE ORDER

#### medipro® RELAX hip disartic pantyhose



Circumference and length (in cm)	Garment size			
	S	M	L	XL
Circumference waist, cT	60 – 70	70 – 85	85 – 100	90 – 115
Circumference thigh, cG	45 – 55	50 – 65	55 – 70	65 – 80

In addition to the garment size we also need:

- The proper seam direction:  
**Seam direction type A** (hip disarticulation) or  
**Seam direction type B** (hemipelvectomy)
- Amputated side: left leg, right leg, both legs

#### medipro® RELAX single-legged pantyhose (max. leg length 20 cm)




Art. No. M.870.902

medipro® care neutral soap, 200 ml

For thorough cleansing of the residual limb and the stocking.

Best if ordered together with the stocking.



Circumference and length (in cm)	Garment size			
	S	M	L	XL
Circumference waist, cT	60 – 70	70 – 85	85 – 100	90 – 115
Circumference hip, cH	80 – 95	85 – 100	95 – 115	110 – 125
Circumference thigh, cG	45 – 55	50 – 65	55 – 70	65 – 80

In addition to the garment size we also need:

- Length of residual limb (IG) in cm (max. leg length 20 cm)
- Amputated side: left leg, right leg, both legs




# medipro® fax order

## Residual limb compression stockings

### medipro® RELAX hip disartic / single-legged pantyhose

**FAX ORDER** to your Distributor

#### medipro® RELAX hip disartic pantyhose



- Number of garments \_\_\_\_\_
- Special detergent Art. No. M.870.902 \_\_\_\_\_ (bottles)  
medipro® care neutral soap, 200 ml

■ Please choose the required pantyhose model:


	S <input type="checkbox"/>	M <input type="checkbox"/>	L <input type="checkbox"/>	XL <input type="checkbox"/>
Circumference waist cT (cm)	60-70	70-85	85-100	90-115
Circumference thigh cG (cm)	45-55	50-65	55-70	65-80

■ Seam direction

- Hipdisarticulation  (seam direction type A)
- Hemipelvectomy  (seam direction type B)

- Amputation: left leg  right leg  both sides

#### medipro® RELAX Single-legged pantyhose



- Number of garments \_\_\_\_\_
- Special detergent Art. No. M.870.902 \_\_\_\_\_ (bottles)  
medipro® care neutral soap, 200 ml

■ Please choose the required pantyhose model:

	S <input type="checkbox"/>	M <input type="checkbox"/>	L <input type="checkbox"/>	XL <input type="checkbox"/>
Circumference waist cT (cm)	60-70	70-85	85-100	90-115
Circumference hip cH (cm)	80-95	85-100	95-115	110-125
Circumference thigh cG (cm)	45-55	50-65	55-70	65-80
Max. leg length (cm)	20	20	20	20

■ Residual limb length (IG)

- Left amputation  \_\_\_\_\_ cm (max. length 20 cm)
- Right amputation  \_\_\_\_\_ cm (max. length 20 cm)

medi. i feel better.

#### Your data

Customer No. \_\_\_\_\_

Commission \_\_\_\_\_

Company \_\_\_\_\_

Contact/Function \_\_\_\_\_

Address \_\_\_\_\_

Postcode/Town \_\_\_\_\_

Telephone \_\_\_\_\_

Fax \_\_\_\_\_

Email \_\_\_\_\_

Date \_\_\_\_\_

Stamp

Signature \_\_\_\_\_

**Please note the delivery period of about 10 working days!**

Head-Office medi • medipro® Prothetik

Medicusstr. 1 • D-95448 Bayreuth • Germany • phone +49 921 912 740 • fax +49 921 912 755  
prothetik@medi.de • www.medi.de • www.stolperstein.com

COPY DRAFT: FAXORDER medipro® RELAX hip disartic / single-legged pantyhose

99.804/09.07/3,5'



4 10263981430836



**medi Austria GmbH**

Steinbockweg 34  
A-6026 Innsbruck  
Tel.: +43-512-57 95 15  
Fax: +43-512-57 95 15 45



**medi Belgium bvba**

Voortstraat 79  
B-3580 Beringen  
Tel.: +32-11-24 25 60  
Fax: +32-11-24 25 64



**medi Bayreuth España SL**

C/Canigó 2 - 6 bajos  
E-08901 Hospitalet de  
Llobregat Barcelona  
Tel.: +34-93-2 60 04 00  
Fax: +34-93-2 60 23 14



**medi France**

ZI Charles de Gaulle  
25, Rue Henri Farman  
F-93297 Trembley en France Cedex  
Tel.: +33-1-48 61 76 10  
Fax: +33-1-49 63 33 05



**medi Italia s.r.l.**

Via Guido Rossa 28  
I-40033 Casalecchio di  
Reno-Bologna  
Tel.: +39-051-6 13 24 84  
Fax: +39-051-6 13 29 56



**Mssrs, medi Bayreuth**

Imp. De Artigos Medicos, Lda.  
Rua Eugénio de Castro Rodrigues, n° 2  
P-1700 Lisboa  
Tel.: +351-21-8 43 71 60  
Fax: +351-21-8 47 08 33



**medi Danmark ApS**

Ved Slotsbanken 14  
DK-6270 Tønder  
Tel.: +45-74-72 01 33  
Fax: +45-74-72 01 63



**RSL Steeper Leeds Manufact.** · Unit 7 Hunslet Trading Estate  
Severn Road · LS101B Leeds LS 10 1BL / England  
Tel.: +44-11-32 72 56 80 · Fax: +44-11-32 76 03 89

mediven®

mediortho®

medipro®



I feel better.

Head-Office medi • medipro® Prothetik

Medicusstr. 1 · D-95448 Bayreuth · Germany · phone +49 921 912 740 · fax +49 921 912 755  
prothetik@medi.de · www.medi.de · www.stolperstein.com