

# REVOLUTION SERIES™



Jeremy Burleson, World-class  
Paralympic Pentathlete



A REVOLUTION IN MOTION




# REVOLUTION SERIES™

## renegade/'reni-gād

*adj* 1: having rejected tradition; 2: unconventional

Break away from the pack and experience a revolution in motion with the truly innovative, patent-pending **renegade**, the newest foot from Freedom Innovations. The **renegade** stands alone as the industry's lightest, most responsive shock foot, driving shock absorption in a revolutionary direction. Beyond a conventional shock system, the **renegade** also redirects the energy to enhance the foot's heel dynamics and propel the amputee through the gait cycle, providing unmatched comfort and performance.




### the revolutionary **renegade** foot

The **renegade**, the most innovative prosthetic foot in the industry today:

- Offers incredibly smooth gait at any speed
- Incorporates a full-length, split-shank, floating heel and composite Z-shock
- Offers unmatched energy return compared to other shock systems
- Achieves unprecedented tibial progression level
- Industry's lightest shock absorbing foot system
- Designed as the ideal prosthetic foot for active amputees

# the renegade:

a revolution in shock absorption



the industry's  
most responsive  
shock foot

As the industry's most responsive shock foot, the renegade's redirected energy propels the amputee through the gait cycle, providing unmatched performance.



# REVOLUTION SERIES™

## THE INDUSTRY'S MOST INNOVATIVE PROSTHETICS FOR ACTIVE AMPUTEES

Freedom Innovations is an advanced-technology prosthetics company, creating a revolution in motion as a global innovator and leader in evolving prosthetics. Merging extensive experience in composite design and prosthetic development with the industry's most innovative technologies, Freedom Innovations creates the industry's highest quality, top performing prosthetic products for amputees worldwide.

### Revolution Series™ Selection Guide

Weight (lbs.)	Impact Level			Weight (kg.)
	Low	Moderate	High	
100-115	1	1	2	44-52
116-130	1	2	3	53-59
131-150	2	3	4	60-68
151-170	3	4	5	69-77
171-195	4	5	6	78-88
196-220	5	6	7	89-100
221-255	6	7	8	101-116
256-285	7	8	9	117-130
286-325	8	9		131-147
326-365	9			148-166

#### ORDERING EXAMPLE - NRS1826R

NRS1	Renegade Foot Category	Size	Left/Right
	8	26cm	R



Products Division, Unit 7, Hunslet Trading Estate, Severn Road, Leeds LS10 1BL, United Kingdom  
Tel: +44 (0)870 240 4133 Fax: +44 (0)870 240 4779  
Email: sales@rslsteeper.co.uk Web: www.rslsteeper.com

## A REVOLUTION IN MOTION