

REVOLUTION SERIES™

The NEW renegade LP. Lowest-profile, lightest shock foot.

The **renegade LP** stands alone as **the industry's lightest, most responsive shock foot**, driving shock absorption to a revolutionary level. Featuring a complete shock system for **clearances under 4"**, the **renegade LP** is ideal for the active amputee with a longer residual limb. Moreover, the **renegade**



LP incorporates a full-length sole plate creating natural inversion and eversion. Composite Z-shock technology with an unprecedented tibial progression level further reinforces our commitment to innovation, value and technology.

the revolutionary **renegade LP**.

Featuring the lightest, most responsive and complete shock system for clearances under 4 inches, the renegade LP is ideal for the active amputee with longer residual limb, driving shock absorption to a revolutionary level.

- Incredibly smooth gait at any speed
- Only shock foot for clearances under 4" (range of 3.7-4.7")
- Full-length sole plate and composite Z-shock technology
- Unmatched energy return compared to other shock systems
- Unprecedented tibial progression level
- Easy to cosmetically finish
- Provides a shapelier ankle



FREEDOM INNOVATION

ONCE AGAIN, THE LEADER IN
INNOVATION, VALUE AND TECHNOLOGY.

Freedom Innovations is an advanced technology prosthetics company, creating a **revolution in motion** as a global innovator and leader in evolving prosthetics.

Merging extensive experience in composite design and prosthetic development with **the industry's most innovative technologies**, **Freedom Innovations** creates the industry's highest quality, top performing prosthetic products for amputees worldwide.

Revolution Series™ Selection Guide

Weight (Ibs.)	Impact Level			Weight (kg.)
	Low	Moderate	High	
100-115	1	1	2	44-52
116-130	1	2	3	53-59
131-150	2	3	4	60-68
151-170	3	4	5	69-77
171-195	4	5	6	78-88
196-220	5	6	7	89-100
221-255	6	7	8	101-116
256-285	7	8	9	117-130
286-325	8	9		131-147
326-365	9			148-166

ORDERING EXAMPLE - NRS3627L

	Renegade LP code	Category	Size	Left/Right
NRS	3	6	27	L



Products Division, Unit 7, Hunslet Trading Estate, Severn Road, Leeds LS10 1BL, United Kingdom

Tel: +44 (0)870 240 4133 Fax: +44 (0)870 240 4779

Email: sales@rslsteeper.co.uk Web: www.rslsteeper.com