

Local details



## Custom Silicone Cosmeses



Silicone Clinic  
Unit 20, Kingsmill Business Park  
Chapel Mill Road  
Kingston upon Thames  
Surrey KT1 3GZ

Email: [silicone.shop@rslsteeper.com](mailto:silicone.shop@rslsteeper.com)  
Tel: (+44) 0870 850 8228  
Fax: (+44) 0870 850 8229

RSLIT252 ISSUE 1 13/05/09

[www.rslsteeper.com](http://www.rslsteeper.com)

**spectroMATCH**





## Custom Silicone from RSLSteeper

Experience the very best in Custom Silicone Cosmeses from RSLSteeper. Each cosmesis is sculpted by hand to meet your requirements and has a remarkably realistic appearance and feel. Skin pigmentation, freckles, veins, hair and even tattoos can be incorporated to ensure an authentic finish.

Unlike many similar products, our fabrication techniques involve colouring the silicone material itself prior to sculpting, as opposed to surface painting. This process allows the layering of different shades of silicone to give a translucent finish to the cosmesis, similar to real human tissue.



### The perfect fusion of technology and art

Utilising Spectromatch, our unique electronic colour matching process, we deliver a near perfect match to your individual skin shade. Our Spectromatch system stores up to 3000 colours and offers a far superior colour matching process than traditional swatches. In addition, numerous colour readings can be completed in a manner of minutes, providing near perfect results regardless of lighting conditions.



### We make:

- individual fingers and toes
- partial feet
- lower limb covers
- partial hands
- upper limb covers
- custom liners
- cosmetic restorations