


# 4Seal TFS Liners

- The 4Seal Liner is designed for use in a socket with a distal expulsion valve and provides excellent suspension by means of its four integral seals.
- The Easy Glide surface and flexible nature of the seals make it simple to roll on and negates the need for any lubricant when donning the socket.
- The Umbrellan fabric contained within the silicone as part of the liner matrix, provides both increased durability and pain-relieving effect from phantom limb pain, phantom limb sensation and idiopathic residual limb pain.



# 4Seal TFS Liners



## Medi 4Seal Comfort Sleeve

These knitted sleeves are available in sizes 26 - 50cm, for use where a proximal reduction in residual limb volume occurs, or to provide increased comfort around the socket brim.

The non slip TPE bands on the distal edge hold the sleeve in place when donning the prosthesis, but should always be positioned above the top seal of the liner.




## Order Information

Size	Part Number	Description
26.5 cm	NM891726	medi 4Seal TFS
28 cm	NM891728	medi 4Seal TFS
30 cm	NM891730	medi 4Seal TFS
32 cm	NM891732	medi 4Seal TFS
34 cm	NM891734	medi 4Seal TFS
36 cm	NM891736	medi 4Seal TFS
38 cm	NM891738	medi 4Seal TFS
40 cm	NM891740	medi 4Seal TFS
45 cm	NM891745	medi 4Seal TFS
50 cm	NM891750	medi 4Seal TFS

Size	Part Number	Description
26 - 32 cm	N1880006	Comfort Sleeve - Size S
34 - 40 cm	N1880007	Comfort Sleeve - Size M
45 - 50 cm	N1880008	Comfort Sleeve - Size L

## Technical Information

Material	Transparent Silicone
Matrix	Continuous fabric within the silicone
Length	40cm
Amputation level	Transfemoral
Items included	Instruction manual, drying stand, Umbrellan pain relief diary and pain scale



Prosthetics Division  
Unit 7  
Hunslet Trading Estate  
Severn Road  
LEEDS LS10 1BL

Tel: +44 (0) 870 240 4133  
Email: customerservices@steepergroup.com

[www.steepergroup.com](http://www.steepergroup.com)

RSLIT353 Issue 2 2015

