

**19**

**TrueFinish™  
Colour Shades**



## Elegance Extra Length Silicone and PVC Cosmetic Gloves

- Extra Length as standard in 7 styles
- Recommended for Steeper, Otto Bock, Motion Control and Hosmer adult and youths, myo electric and mechanical hand ranges
- A choice between the Standard colour range or the advanced micropigmentation TrueFinish™ for added realism and depth

Available in 19 colour shades, 5 sizes and 7 styles

## Elegance Silicone

By using Tufsil 3™ our innovative polymer, Elegance Silicone gloves offer the toughness associated with PVC but provide the ease of cleaning provided by silicones.

In addition, Tufsil 3™ allows additional stretch, reducing the energy expended by the hand motor each time it opens. Properties of this polymer include improved tear resistance proven by clinical evaluation and up to 15% stretch, to allow even bulbous sockets to be accommodated.

All Elegance Silicone gloves are easy to fit and remove, allowing quick access to the hardware for maintenance and adjustment. The Silicone gloves are simple to keep clean using soap and water.

## Elegance PVC

Elegance PVC gloves provide a cost effective solution for electric and mechanical hands when the extra toughness of PVC is required.



Size	Approx Length	Application	Left Silicone	Right Silicone	Left PVC	Right PVC	Colour Ranges
<b>Standard Colours</b>							
6 ¾	460mm	Youths	SG310	SG311	CG310	CG311	E0 - E18
7	520mm	Female / Youths	SG308	SG309	CG308	CG309	E0 - E18
7 ¼	645mm	Youths	SG306	SG307	CG306	CG307	E0 - E18
7 ¼	645mm	Female	SG304	SG305	CG304	CG305	E0 - E18
7 ¾	620mm	Female	SG312	SG313	CG312	CG313	E0 - E18
7 ¾	645mm	Male	SG302	SG303	CG302	CG303	E0 - E18
8 ¼	645mm	Male	SG300	SG301	CG300	CG301	E0 - E18
<b>TrueFinish™ Colours</b>							
6 ¾	460mm	Youths	TSG310	TSG311	TG310	TG311	E0 - E18
7	520mm	Female / Youths	TSG308	TSG309	TG308	TG309	E0 - E18
7 ¼	645mm	Youths	TSG306	TSG307	TG306	TG307	E0 - E18
7 ¼	645mm	Female	TSG304	TSG305	TG304	TG305	E0 - E18
7 ¾	620mm	Female	TSG312	TSG313	TG312	TG313	E0 - E18
7 ¾	645mm	Male	TSG302	TSG303	TG302	TG303	E0 - E18
8 ¼	645mm	Male	TSG300	TSG301	TG300	TG301	E0 - E18

## 19 TrueFinish™ Colour Shades



### Ordering Examples

TSG300/E5

Orders a 8 ¼ Left Silicone Glove in TrueFinish colour E5

CG300/E5

Orders a 8 ¼ Left PVC Glove in Standard Finish colour E5

Swatches are available on request:

E24654 E0 - E18 Standard colours

E24655 E0 - E18 TrueFinish colours