



19
TrueFinish™
Colour Shades



Skinergy Plus Leg Covers

Lifelike, durable and easy to fit, Skinergy Plus offers a fantastic, off-the-shelf, cosmetic finish for the below knee amputee.

Available in 19 skin shades and featuring a NEW silicone formulation for a more natural finish.





Features:

- Integrated reinforced mesh
- 19 colour shades
- Seamless design
- Easy to clean
- Available in closed or split toe
- Adult foot sizes 22-30 cm
- Two calf sizes

Split Toe Part Number	Closed Toe Part Number	Foot Size (cm)	Ankle Circ.		Mid Calf Circ.		Top of Cover Circ.		Cover Length From Floor to -	
			Min	Max	Min	Max	Min	Max	Mid Calf	Cover Top
TLC22LS/E*	TNS22LS/E*	22-23	20	25	28	35	26	32.5	32	55
TLC22RS/E*	TNS22RS/E*	22-23	20	25	28	35	26	32.5	32	55
TLC22LL/E*	TNS22LL/E*	22-23	20	25	32	40	30	37.5	32	55
TLC22RL/E*	TNS22RL/E*	22-23	20	25	32	40	30	37.5	32	55
TLC24LS/E*	TNS24LS/E*	24-25	22	27.5	31	39	29	36	34.5	57
TLC24RS/E*	TNS24RS/E*	24-25	22	27.5	31	39	29	36	34.5	57
TLC24LL/E*	TNS24LL/E*	24-25	22	27.5	35	44	33	41	34.5	57
TLC24RL/E*	TNS24RL/E*	24-25	22	27.5	35	44	33	41	34.5	57
TLC26LS/E*	TNS26LS/E*	26-27	22	28	32	40	30	37.5	37	61
TLC26RS/E*	TNS26RS/E*	26-27	22	28	32	40	30	37.5	37	61
TLC26LL/E*	TNS26LL/E*	26-27	22	28	37	46	35	44	37	61
TLC26RL/E*	TNS26RL/E*	26-27	22	28	37	46	35	44	37	61
TLC28LS/E*	TNS28LS/E*	28-30	24	30	34	42.5	32	40	39.5	66
TLC28RS/E*	TNS28RS/E*	28-30	24	30	34	42.5	32	40	39.5	66
TLC28LL/E*	TNS28LL/E*	28-30	24	30	39	49	37	46	39.5	66
TLC28RL/E*	TNS28RL/E*	28-30	24	30	29	49	37	46	39.5	66

19 TrueFinish™ Colour Shades



Part Number	Description
NCT9740	Skinergy Plus TrueFinish 25mm Punchtool
E24655	Skinergy Plus TrueFinish Swatch Kits
E26625/E*	Skinergy Plus TrueFinish Patch Kits
E26626/E*	Skinergy Plus TrueFinish Spare Parts

Note: E* denotes colour. For example, TNS22LS/E5

Prosthetics Division , Unit 7,
Hunslet Trading Estate, Severn Road
Leeds, LS10 1BL
Tel: +44 (0) 870 240 4133
Email: customerservices@steepergroup.com
www.steepergroup.com