



SENATOR™ PROSTHETIC FOOT

SOPHISTICATED
SIMPLICITY.



The Senator Prosthetic Foot is a lightweight, economical, energy-returning product that shares composite engineering features with our advanced Freedom Series™ prosthetic feet. It is designed to meet the needs of low to moderately active amputees seeking a prosthetic foot that will allow them to ambulate at varied speeds and participate in daily activities such as bicycling, golfing, hiking, or even jogging.



AN EXCELLENT VALUE AT
AN AFFORDABLE PRICE

The Senator Prosthetic Foot combines sophisticated carbon fibre with an integral connector for a lightweight, durable, dynamic-response foot that is simple to fit and finish. It is sophisticated simplicity at its best, delivering substantial energy return and smooth rollover.

MAJOR BENEFITS

- Lightweight, energy-return design facilitates smooth, easy stride so the user can walk further, faster using less energy
- Integral male pyramid attachment simplifies the fitting process
- Simulated ankle motion provides a stable base of support throughout stance phase
- Foot is not pre-bonded to the foot shell, allowing access to the foot module and simplifying the use of heel wedges to improve alignment
- Six stiffness categories and low or moderate impact levels enhance customisation to ensure ultimate comfort and stability
- Ideal for low to moderately active amputees weighing up to 300 lbs (136 kg), enhancing the range of amputees who can utilise this technology
- Low maintenance, durable design provides peace of mind to the clinician and the user



SENATOR™ PROSTHETIC FOOT



CATEGORY SELECTION CHART

Weight (lbs)	Impact Level		Weight (kg)
	Low	Medium	
100-120	1	2	45-55
121-145	2	3	56-66
146-175	3	4	67-80
176-210	4	5	81-96
211-250	5	6	97-114
251-300	6		115-136

ORDERING EXAMPLE - NVS1627L

	Foot Type	Category	Size	Left/Right
NVS	1	6	27	L

PRODUCT CLEARANCE & WEIGHT

Size (cm)	Clearance	Weight	
		Foot	Shell
22-25	4.5 in	362g	127g
26-28	5.0 in	390g	159g
29-30	5.4 in	418g	219g



SPECIFICATIONS

Size offering - 22-30cm
Carbon fibre foot warranty - 24 months
Foot shell warranty - 6 months

Low Impact - Average to brisk walking, golf, occasional hiking on uneven trails, or possibly light jogging

Medium Impact - More active, such as fast jogging, carrying heavy loads; engaging in tennis, golf, light jogging or hiking on uneven surfaces on regular basis

RSL STEEPER
Products Division
Unit 7, Hunslet Trading Estate
Severn Road, Leeds LS10 1BL
United Kingdom
Tel: +44 (0)870 240 4133
Fax: +44 (0)870 240 4779
Email: sales@rslsteeper.co.uk
Web: www.rslsteeper.com