

SKINERGY PLUS™

Prosthetics Division



Skinerger Plus leg covers

More durable and realistic than ever



 **RSLSTEEPER**
Listening to you, working with you

Skinergy Plus from



The natural choice

Skinergy leg covers are a popular choice the world over – their realistic appearance and convenience enhance the everyday lives of thousands.

Now, Skinergy Plus takes this experience to a new level.

For people who like to lead outdoor and demanding lifestyles, the advanced materials and increased durability of Skinergy Plus will help them to get more out of life than ever.

Easy to fit, enhanced toughness and improved realism make Skinergy Plus our most impressive leg cover yet – the natural recommendation for your patients.




Increased durability, more shade options

Skinergy Plus sets new standards for silicone cosmeses:

- **Integrated reinforcing mesh** and **advanced silicone material** for a more durable finish
- Now available in **19 colour shades** with unique **True Finish™** micro-pigmentation for authentic realism

Skinergy Plus still offers the great features that Skinergy customers love:

- Seamless design
- Easy to clean
- Available in closed or split foot
- Adult foot sizes 22cm to 30cm
- Two calf sizes



Part No. Split Toe	Part No. Closed Toe	Foot Size	Ankle Circ.		Mid Calf Circ.		Top of Cover Circ.		Cover Length from Floor to - Mid Calf Cover Top	
			Min.	Max.	Min.	Max.	Min.	Max.		
TLC22LS/E*	TNS22LS/E*	220 / 230	20	25	28	35	26	32.5	32	55
TLC22RS/E*	TNS22RS/E*	220 / 230	20	25	28	35	26	32.5	32	55
TLC22LL/E*	TNS22LL/E*	220 / 230	20	25	32	40	30	37.5	32	55
TLC22RL/E*	TNS22RL/E*	220 / 230	20	25	32	40	30	37.5	32	55
TLC24LS/E*	TNS24LS/E*	240 / 250	22	27.5	31	39	29	36	34.5	57
TLC24RS/E*	TNS24RS/E*	240 / 250	22	27.5	31	39	29	36	34.5	57
TLC24LL/E*	TNS24LL/E*	240 / 250	22	27.5	35	44	33	41	34.5	57
TLC24RL/E*	TNS24RL/E*	240 / 250	22	27.5	35	44	33	41	34.5	57
TLC26LS/E*	TNS26LS/E*	260 / 270	22	28	32	40	30	37.5	37	61
TLC26RS/E*	TNS26RS/E*	260 / 270	22	28	32	40	30	37.5	37	61
TLC26LL/E*	TNS26LL/E*	260 / 270	22	28	37	46	35	44	37	61
TLC26RL/E*	TNS26RL/E*	260 / 270	22	28	37	46	35	44	37	61
TLC28LS/E*	TNS28LS/E*	280 / 300	24	30	34	42.5	32	40	39.5	66
TLC28RS/E*	TNS28RS/E*	280 / 300	24	30	34	42.5	32	40	39.5	66
TLC28LL/E*	TNS28LL/E*	280 / 300	24	30	39	49	37	46	39.5	66
TLC28RL/E*	TNS28RL/E*	280 / 300	24	30	39	49	37	46	39.5	66

How to order

You can order Skinergy Plus leg covers using the above reference numbers.

For more information or to place an order, please contact RSLSTEEPER:

Telephone: 0870 240 4133
 Email: sales@rslsteeper.co.uk
www.rslsteeper.com

Description	Part No.
Skinergy Plus True Finish 25mm Punchtool	NCT9740
Skinergy Plus True Finish Swatch Kits	E24658
Skinergy Plus True Finish Patch Kits	E26625/E*
Skinergy Plus True Finish Spare Kits	E26626/E*

E* denotes colour. For example, TNS22LS/E5