

# Plié 3 MPK Package Order Form

## Included in the package:

- Plié 3 knee, including 5 year warranty & the following items:
  - Wireless USB adapter
  - Lithium-ion battery case
  - Software CD
  - Lithium-ion battery x2
  - Battery charger with car adapter
  - Air pump with hose adapter
  - Smooth hose adapter for air pump
- Preferred designer fairing
- Choice of foot - Kinterra, Sierra or Highlander, Maverick Xtreme, Maverick Xtreme AT and Maverick Comfort AT
- Set of adapters

## How to order

Step 1	Step 2	Step 3	Step 4 - Successful Trial Only
Submit a PO using part number: <b>KIT-PLIEMPK</b> and your completed order form.	Once the PO and patient form has been received, this will generate an order for a 30 day trial period for the <b>KIT-PLIEMPK</b> .  You will not be charged or invoiced at this stage.	Following a successful trial, call Customer Services on 0870 240 4133 to obtain an RTA number and return the demo Plié 3 knee (including accessories and case) only.  If the trial was unsuccessful, call Customer Services on 0870 240 4133 to obtain an RTA number and return the full kit.	After contacting Customer Services your order will then be finalised. The definitive Plié 3 knee and the preferred designer fairing stated on this form will be delivered and invoiced.

## Patient Information:

Name:	Age:	Gender:	M	F	
Weight (kg) :	Height:	Skin tone:	Caucasian:	Mid Brown:	Brown:
Amputation Level (TF, TT, KD, HD...):	Left:	Right:	Size (22-31cm):		
Impact level (see descriptions overleaf):	Low:	Moderate:	High:		
Notes:					
Prosthetist:	Branch / Contractor:				


## Plié 3 Knee

Pyramid Adapter:  Threaded Adapter:  *Please note that only a pyramid adapter is supplied with trial knees.*

## Designer Fairings

Please select the preferred designer fairing. This will be supplied on the definitive Plié 3 knee upon completion of a successful evaluation and will not be shipped with the demonstration unit. The demonstration unit will be supplied in Standard Black.

- Standard Black
- Proud Marines
- Captivating Copper
- Radical Red
- The Patriot
- Pink Passion
- Forest Camo
- Fabulous Flames
- Tech Metal
- Blazing Blue
- Pearl White



## Adapters

The following adapters are included within the MPK package. For alternative adapters, please place a separate order.

### Set 1

Description	Part Number	Weight Rating	Material
Tube Clamp, 30 mm, Female Pyramid	NACC-00-12500-00	136kg	Titanium
Tube w/Titanium Female Pyramid, 30 mm x 400 mm	NACC-00-12800-00	136kg	Titanium
Pyramid, 4 Hole, Female	NACC-00-13200-00	136kg	Titanium

### Set 2

Description	Part Number	Weight Rating	Material
Tube Clamp, 30 mm, Female Pyramid	NACC-00-12500-00	136kg	Titanium
Tube w/Titanium Female Pyramid, 30 mm x 400 mm	NACC-00-12800-00	136kg	Titanium
Base, 3 Prong, Rotating Pyramid	NACC-00-13800-00	136kg	Stainless Steel
Pyramid, Rotating, Female	NACC-00-14000-00	136kg	Titanium

## Foot Selection



Category Selection - Kinterra			Category Selection - Sierra, Highlander, Maverick Xtreme, Maverick Xtreme AT and Maverick Comfort AT			
Weight (kg)	Impact Level		Weight (kg)	Impact Level		
	LOW	MOD		LOW	MOD	HIGH
44-52	1	1	44-52	1	1	2
53-59	1	2	53-59	1	2	3
60-68	2	3	60-68	2	3	4
69-77	3	4	69-77	3	4	5
78-88	4	5	78-88	4	5	6
89-100	5	6	89-100	5	6	7
101-116	6	7	101-116	6	7	8
117-125	7		117-130	7	8	9
			131-147	8	9	
			148-166	9		

## Side

Left: Right:

## Toe Configuration<sup>2</sup>

Regular Toe: Sandal Toe:

## Size (cm)

22 23 24 25 26 27 28 29 30 31

## Stiffness Category

1 2 3 4 5 6 7 8<sup>1</sup> 9<sup>1</sup>

<sup>1</sup>Sierra, Highlander, Maverick Xtreme, Maverick Xtreme AT & Maverick Comfort AT only.

<sup>2</sup>Sandal toe available in sizes 22-28 for Kinterra, Sierra & Highlander. Sandal toe available in sizes 23-28 for all Maverick feet.

### Impact Level Descriptions

**Low:** Daily activities include mostly ground level walking, moving around in the home and the community.

**Moderate:** Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

**High:** Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

## Plié 3 Warranty

5 years

The extended warranty on the Plié 3 from 3 to 5 years is automatically included within the MPK package. Please note, the Plié 3 knee must be returned for annual servicing to maintain the warranty. A loaner unit will be provided during this time.

For warranty terms on all other items, please refer to the Extended Warranty - KIT-PLIEMPK information.

Please return your completed order form, PO and reference to KIT-PLIEMPK to [customerservices@steepergroup.com](mailto:customerservices@steepergroup.com)

Tel: +44 (0) 870 240 4133 Fax: +44 (0) 870 240 4779

[www.steepergroup.com](http://www.steepergroup.com)

STPPR031 Issue 4 2018